

Supplementary Table IV. Primer sequence used for q-RT-PCR validation. Primer sequence was designed for mature miRNA sequence and universal reverse primer was used

Sr. no.	Mir-RNA (forward)	Sequence
1	hsa-miR-193b-5p	CGGGUUUUGAGGGCGAGAUGA
2	hsa-miR-326	CCUCUGGGCCCUUCCUCCAG
4	hsa-miR-1307-5p	UCGACCGGACCUCGACCGGCU
5	hsa-miR-671-5p	AGGAAGCCCUGGAGGGGCUGGAG
7	Reverse primer	Poly T tail universal primer